

Cold Mounting



Description

Cold Mounting system offer excellent pore impregnation in a compact format. As a result, edge retention is enhanced and friable samples are supported during grinding and polishing.

Features

- **High Efficiency :** The rugged vacuum pump supplies ample vacuum pull down to quickly evacuate trapped air from any porous specimen while the patented, high strength plastic chamber maintains vacuum throughout the impregnation period.
- **Optimized Pore Impregnation :** The unique dispensing mechanism with a built in rotating turn table. Allows pouring epoxy while under vacuum.